

MODULAR ARTICULATOR SYSTEM



1210 AS Basic Module

Includes:

- 1701 AR PCH Articulator (1.5mm analogs)
- 2094 AA Articulator Support Legs (set)
- 4000 FB **Pana-Mount** Face-bow
- 4070 FB **Bio-Esthetic** Level Gauge
- 8210 AA Molded Plastic Case
- 2280 ME Mandibular Mounting Stand
- 2850 MP Magna-Split II System
- 7172 DP **Bite-Tray** Registration Plates (20)
- 4405 DP **Bite-Tab** Compound Discs (180 Disks)
- 5401 EM Basic Instruction Manual
- 9850 MP Pana-Lube Silicone Spray

1211 AS Basic Module

Includes:

- 1801 AR PSH Articulator (1.5mm analogs)
- 2094 AA Articulator Support Legs (set)
- 4000 FB **Pana-Mount** Face-bow
- 4070 FB **Bio-Esthetic** Level Gauge
- 8210 AA Molded Plastic Case
- 2280 ME Mandibular Mounting Stand
- 2850 MP Magna-Split II System
- 7172 DP **Bite-Tray** Registration Plates (20)
- 4405 DP **Bite-Tab** Compound Discs (180 Disks)
- 5401 EM Basic Instruction Manual
- 9850 MP Pana-Lube Silicone Spray

The basic module includes all of the essential items needed to mount casts in centric relation and set the protrusive condylar analog guide with a unique intraoral "check bite" jaw positional record. The molded plastic case holds the articulator, face-bow and accessories plus room for mounted casts. The basic module routinely comes with 1.5mm Bennett analogs. The module comes with magnetic mounting plate system.

MODULAR ARTICULATOR SYSTEM



1230 AS Advanced Module

Includes:

- 1701 AR PCH Articulator (1.5mm analogs)
- 2094 AA Articulator Support Legs (set)
- 4000 FB **Pana-Mount** Face-bow
- 4070 FB **Bio-Esthetic** Level Gauge
- 8210 AA Molded Plastic Case
- 2280 ME Mandibular Mounting Stand
- 7172 DP **Bite-Tray** Registration Plates (20)
- 4405 DP **Bite-Tab** Compound Discs (180 Disks)
- 7600 AP Axi-Path III Recorder (Macro-Micro adjust Axis Arms)
- 2250 ME Axis Mounting System
- 2804 MP Magna-Split Mounting Plate (for AMS)
- 2850 MP Magna-Split II System
- 2105 AA 0.5mm Motion Analogs
- 2110 AA 1.0mm Motion Analogs
- 2120 AA 2.0mm Motion Analogs
- 2125 AA 2.5mm Motion Analogs
- 5401 EM Basic Instruction Manual
- 5403 EM Advanced Instruction Manual
- 9850 MP Pana-Lube Silicone Spray

1231 AS Advanced Module

Includes:

- 1801 AR PSH Articulator (1.5mm analogs)
- 2094 AA Articulator Support Legs (set)
- 4000 FB **Pana-Mount** Face-bow
- 4070 FB **Bio-Esthetic** Level Gauge
- 8210 AA Molded Plastic Case
- 2280 ME Mandibular Mounting Stand
- 7172 DP **Bite-Tray** Registration Plates (20)
- 4405 DP **Bite-Tab** Compound Discs (180 Disks)
- 7600 AP Axi-Path III Recorder (Macro adjust Axis Arms)
- 2250 ME Axis Mounting System
- 2804 MP Magna-Split Mounting Plate (for AMS)
- 2850 MP Magna-Split II System
- 2105 AA 0.5mm Motion Analogs
- 2110 AA 1.0mm Motion Analogs
- 2120 AA 2.0mm Motion Analogs
- 2125 AA 2.5mm Motion Analogs
- 5401 EM Basic Instruction Manual
- 5403 EM Advanced Instruction Manual
- 9850 MP Pana-Lube Silicone Spray

The Advanced Module includes all items in the Basic Module, as well as the Axi-Path III Recorder and Axis Mounting System. The Axi-Path III Recorder provides the means to quickly and accurately record a patient's mandibular movements on the transverse (hinge) axis. The Axis Mounting System is used for true hinge axis transfer mounting as well as arbitrary face-bow mountings. Also included is the full set of motion analogs (0.5, 1.0, 1.5, 2.0, 2.5mm).