

MOUNTING PLATES

2850 MP *Magna-Split II™ System*

Includes:

2 Split-Cast Magnetic Mounting Plates, 20 Plastic Index Plates



2860 MP *Magna-Split™ Index Plates* (20)

2861 MP *Magna-Split™ Index Plates* (50)


2862 MP *Magna-Split™ Index Plates* (100)

These pre-formed plastic Index Plates attach directly to the magna-Split Mounting Plate for quick and easy mounting of study casts.



2804 MP *Magna-Split™ Mounting Plate*

The mounting plate comes with removable rare earth magnet and two stainless steel cap screws for fixing the plate to upper or lower frame of the articulator. (Picture shows magnet in place.)



2706 MP *Plastic Mounting Plates* (10)

2707 MP *Plastic Mounting Plates* (20)

2708 MP *Plastic Mounting Plates* (50)

These Plastic Mounting Plates are inexpensive standard mounting plates. They can be used with several brands of articulators including Panadent, Whip-Mix, TMJ and Stuart.

